





HQS product range X-Cath



X-Cath device has been specifically designed to adapt HQS Introducer into multiple accesses. X-Cath device ensures a sealing all over the procedure, using simultaneously 2,3 or 4 tools with the same introducer.

References	2064-X-Cath-2	2064-X-Cath-3	2064-X-Cath-4
Inner view			
Packaging	Unitary Sterile		

References : 2064-X-CATH-A with A the number of accesses (2, 3 or 4)

- Use exclusively with HQS Aalseal introducers 2064-HQS.
- Single use,
- FDA cleared, **CE**0086 marked.



Innovation Show case
2010 Charing Cross 2011



ALSEAL
Temis Santé - Bioparc
2 rue Paul Milleret - 25000 Besançon - France

Phone: +33 (0)3 81 61 69 93
aalseal@alcis.net

Fax: +33 (0)3 81 53 47 65
www.alcis.net/aalseal

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. Always refer to the appropriate instructions for use for complete clinical instructions. Non contractual document. Aalseal reserves the right, without prior notice, to modify the products in order to improve their quality. Always refer to the Instruction For Use and labels of the product before use. WARNING: Applicable laws restrict these products to sale by or on the order of a physician. Legal manufacturer : Alcis. Notified body 0086 for CE mark : BSI. Unless specifically identified as "NOT CE MARKED AND NOT FDA CLEARED.", all the references numbers mentioned on this document are FDA cleared and CE marked according to European council directive 93/42/EEC on medical devices and its relatives.